

MATERIAL SAFETY DATA SHEET

FISH Wash Buffer 2

Pub No.: M224

Rev.: 10/04/2012

SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

PRODUCT NAME: FISH Wash Buffer 2
PART NO: 2010-00-2
MANUFACTURER: SciGene Corporation
 470F Lakeside Drive
 Sunnyvale, CA 94085 USA
EMERGENCY PHONE: Poison Control
 In U.S. 800-222-1222
 Outside U.S. Contact your local poison control center

SECTION 2: HAZARDS IDENTIFICATION

Emergency Overview

OSHA Hazards: Irritant

GHS Classification: Acute toxicity (Category 5)
 Serious eye damage (Category 1)

GHS Label elements, including precautionary statements

Pictogram:



Signal word: Danger

Hazard statements: H303 May be harmful if swallowed
 H318 Causes serious eye damage

HMIS Classification

Health hazard: 2

Flammability: 0

Physical hazards: 0

NFPA Rating

Health hazard: 2

Fire: 0

Reactivity Hazard: 0

Potential Health Effects

Inhalation: May be harmful if inhaled. Causes respiratory tract irritation.

Skin: May be harmful if absorbed through skin. Causes skin irritation.

Eyes: Causes eye irritation.

Ingestion: May be harmful if swallowed.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

NAME	CAS NO	%
α -[(1,1,3,3-Tetramethylbutyl)phenyl]- ω -hydroxy-poly(oxy-1,2-ethanediyl)	9036-19-5	0.1
2.0 X Saline Sodium Citrate		99.7%

There are no ingredients or additional ingredients present which, within the current knowledge of the supplier and in the concentrations applicable, are classified as hazardous to health or the environment and hence require reporting in this section.

SECTION 4: FIRST AID MEASURES

General advice: Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.

If inhaled: If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact: Wash off with soap and plenty of water. Consult a physician.

In case of eye contact: Flush eyes with water as a precaution.

If swallowed: Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician

MATERIAL SAFETY DATA SHEET

FISH Wash Buffer 2

Pub No.: M224

Rev.: 10/04/2012

SECTION 5: FIRE-FIGHTING MEASURES

Flammability:	Not flammable or combustible
Products of combustion:	Carbon oxides
Suitable Extinguishing Media:	Regular dry chemical, carbon dioxide, water spray, regular foam
Protective Equipment:	Wear self contained breathing apparatus (SCBA) for fire fighting if necessary.
Sensitivity to Mechanical Impact:	Not sensitive
Sensitivity to Static Discharge:	No

SECTION 6: ACCIDENTAL RELEASE MEASURES

Personal precautions:	Use personal protective equipment. Avoid breathing vapors or mist. Ensure adequate ventilation. Evacuate personnel to safe areas.
Environmental precautions:	Do not let product enter drains.
Methods and materials for containment and clean up:	Soak up with inert absorbent material. Collect spilled material in appropriate container for disposal.

SECTION 7: HANDLING AND STORAGE

Handling:	Avoid breathing vapors. Wash thoroughly after handling.
Storage:	Store in accordance with local regulations. Keep container closed.

SECTION 8: EXPOSURE CONTROLS / PERSONAL PROTECTION

Hygiene Measures:	Handle in accordance with good industrial hygiene and safety practices. Wash hands, forearms, and face thoroughly after handling chemical products.
Skin and body protection:	Personal protective equipment for the skin and body should be selected based on the task being performed and should be approved by a specialist before handling this product.
Eyes:	Goggles or safety glasses with eyeshields approved to an appropriate government standard.
Hands:	Handle with chemical resistant protective gloves.
Respiratory:	Emmissions from ventilation or work process equipment should be checked to ensure they comply with the requirements of environmental protection legislation. In some cases, the use of fume scrubbers, filters, or properly fitted, air-purifying or air-fed respirators may be necessary.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

Appearance:	Clear liquid
Color:	Colorless
Physical state:	Liquid
pH:	Neutral
Boiling/condensation point:	N/A
Melting /freezing point:	N/A
Solubility:	Soluble in water
Odor:	none
Viscosity:	N/A
Explosion limits:	N/A
Flash point:	N/A
Density:	N/A
Flammability:	Non-flammable

SECTION 10: STABILITY AND REACTIVITY

Stability and reactivity:	Stable under normal conditions.
Conditions to avoid:	No data available.
Incompatibility with various substances:	Strong oxidizing agents.
Hazardous decomposition products:	Carbon monoxide, carbon dioxide.
Hazardous Polymerization:	Will not occur under normal conditions.

SECTION 11: TOXICOLOGICAL INFORMATION

LD50 not established for this Product	
Acute toxicity:	no data available
Other information on acute toxicity:	
Skin corrosion/irritation:	no data available
Serious eye damage/eye irritation:	no data available.
Respiratory or skin sensitization:	no data available
Germ cell mutagenicity:	no data available

MATERIAL SAFETY DATA SHEET

FISH Wash Buffer 2

Pub No.: M224

Rev.: 10/04/2012

SECTION 11: TOXICOLOGICAL INFORMATION (Continued)

Carcinogenicity	This product is or contains a component that is not classifiable as to its carcinogenicity based on its IARC, ACGIH, NTP, or EPA classification.
IARC:	No component of this product present at levels greater than or equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.
ACGIH	No component of this product present at levels greater than or equal to 0.1% is identified as a Carcinogen or potential carcinogen by ACGIH.
NTP:	No component of this product present at levels greater than or equal to 0.1% is identified as a known or anticipated carcinogen by NTP.
OSHA:	No component of this product present at levels greater than or equal to 0.1% is identified as a carcinogen or potential carcinogen by OSHA.
Reproductive toxicity	no data available
Teratogenicity	no data available
Specific target organ toxicity - single exposure (Globally Harmonized System):	no data available
Specific target organ toxicity - repeated exposure (Globally Harmonized System):	no data available
Aspiration hazard:	May be fatal if swallowed and enters airways.
Potential health effects:	
Inhalation:	May be harmful if inhaled. Causes respiratory tract irritation.
Ingestion:	May be harmful if swallowed.
Skin:	May be harmful if absorbed through skin. Causes skin irritation.
Eyes:	Causes eye irritation.
Signs and Symptoms of Exposure:	Ingestion of large amounts may cause nausea; diarrhea
Synergistic effects:	no data available
Additional Information:	RTECS: Not available

SECTION 12: ECOLOGICAL INFORMATION

Persistence and degradability:	no data available
Bioaccumulative potential:	no data available
Mobility in soil:	no data available
PBT and vPvB assessment:	no data available
Other adverse effects:	no data available
Toxicity:	no data available

SECTION 13: DISPOSAL CONSIDERATIONS

Waste disposal: Do not pour down drains. Dispose of material in accordance with all Federal, State and Local regulations.

Refer to **SECTION 7: HANDLING AND STORAGE** and **SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION** for additional handling information and protection of employees.

SECTION 14: TRANSPORT INFORMATION

US DOT Information (49 CFR):	Not dangerous goods
TDG Information:	Not dangerous goods
IATA Information:	Not dangerous goods

SECTION 15: REGULATORY INFORMATION

OSHA Hazards:	Irritant
SARA 302 Components:	SARA 302: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.
SARA 313 Components:	SARA 313: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 313.
SARA 311/312 Hazards:	Acute health hazard
California Prop. 65 Components:	This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

SECTION 16: OTHER INFORMATION

Notice to reader

DISCLAIMER: This Material Safety Data Sheet is offered without charge to the clients of SciGene. Data is the most current available to SciGene at the time of preparation and is issued as a matter of information only, no warranty as to its accuracy or completeness is expressed or implied.